Orlando Bariatric

Chairs, engineered for care.









Born out of human understanding and clinical know-how, our therapeutic chairs are engineered and built to support the body where it most needs support. Adapting to the individuals' needs and allowing weight distribution on pressure points to be spread; our chairs alleviate the stress on the user and the strain on the carer.

Contents

- 3 Introduction
- 4 Overview
- 6 Features
- 8 4 Principles of Effective Clinical Seating
- 9 Infection Control
- 10 Parts & Accessories
- 11 Technical Specifications




Orlando Bariatric

The Orlando™ Bariatric Riser Recliner has been designed to combine the comfort and usability associated with regular bariatric recliners with the clinical seating principles, quality, and robustness that users and caregivers have come to expect from Seating Matters chairs.

The combination of comfort, usability, robustness, and clinical principles offers a unique, clinically advanced bariatric recliner that is equally at home in a private or acute setting. Orlando offers clients the ability to relax and feel comfortable both mentally and physically when in the chair.

Orlando Bariatric has enhanced mobility features to enable it to meet the needs of longer-term care and acute facilities. Easier to manoeuvre, it reduces strain on the caregiver.





04

Orlando Bariatric Features

4

Forward Tilt & Rise +

30° rise angle to suppor independent transfers.4



GoFlat™ Function

Removable Back

Waterfall back comes as standard for optimal postural support and comfort. The entire back is fully removable for ease of delivery and transport.



High Weight Capacity

40-70st*/ 560lb-980lb*/ 254-445kg* weight capacit 20-22st / 280-308lb / 127-140kg leg lift capacity.



Tilt in Space

Tilt in Space with indpendent rise, recline and leg elevation for comfort, repositioning and pressure loading¹ in line with clinical guidelines.²



7

Chair Coverings



Dartex® (As standard)

Waterproof, breathable & + multi-stretch material used for skin contact areas.

Excellent for pressure care.³



Vinyl (As standard)

Waterproof covering which i hard wearing making it ideal for a multi-user environment



Wipeable Fabric Option

fabric which is treated with a Teflon stain resistant finish for easier cleaning.



- 9 Simple Intuitive Handset+
- 10 Pressure Relieving Envelo Cushion



06

Stand Aid / Hoist Compatibility

Simple, safe and effective






4 Principles of Effective, Clinical Seating



Load the Body

Ensure full body is in contact with the chair to redistribute pressure from bony prominences.



2

Provide Postural Support

Accommodate fixed deviations and correct flexible deviations to prevent further progression of deformities and avoid uneven pressure redistribution.



Effective Repositioning

Guidelines recommend that seated individuals are repositioned frequently to increase blood flow to skin. Functions such as tilt in space and back angle recline make this easier.





Appropriate Surface

All weight-bearing parts of the chair should facilitate immersion, envelopment and support of the patient's body. The surfaces of Seating Matters chairs are made of high performance pressure distributing materials with a multi-stretch Dartex covering.



Infection Control

Infection prevention is vitally important in the home, hospital and healthcare environments. Bacteria and germs often contribute to the development and spread of infections. They can potentially be very harmful so this is important to control and prevent, especially in areas such as hospitals, where patients may be vulnerable and susceptible to infection.



At Seating Matters we take infection control very seriously. Our chairs have been specifically designed to decrease entrapment areas where bacteria can survive. The chairs can generally be cleaned using the same process and materials as hospital grade beds and mattresses.

Dartex is a proven multi-stretch patient contact surface with bio-static (anti-mycotic) coatings. It is waterproof, stain resistant, flame retardant and has multi-directional stretch for pressure redistribution. It is used on the seat cushion, calf panels, lateral supports and seat back.

08

Upgrades & Accessories





Build quality

It takes continuous, careful planning to ensure our chairs are functional as well as being aesthetically pleasing. Our design team are adamant that the chairs need to be designed in such a way that they fit into any environment. Each and every part of a Seating Matters chair has been engineered and devised to offer optimum quality and support for the user.

Ease of use

10


All Seating Matters chairs can be adjusted without tools in less than one minute.
Useful to easily and efficiently adjust the chair to suit the user or to readjust for use within a multi-user environment.


References


 Staas, W.E., Cioschi, H.M. (1991). Pressure sores: a multifaceted approach to prevention and treatment. Western Journal of Medicine, 154, 539-544.

Laterals

- 2. Tierney, M. (2022). The Clinician's Seating Handbook, A Reference Guide for Clinical Seating Provision, 6, 62-63
- Trelleborg (2019). Retrieved May 7, 2022. Are You Sitting Comfortably? https://www.trelleborg.com/engineeredcoated-fabrics/-/media/engineered-fabrics/dartex/dartexpdfs/2019-are-you-sitting-comfortably-ecf-dartex.pdf
- National Institute for Health and Care Excellence (NICE) (2009). Retrieved May 10, 2023. Rehabilitation after critical illness in adults [CG83]. https://www.nice.org.uk/guidance/ cn83







Technical Specifications

Recline Leg Rest Eleva	tion Rise Tilt In Space		
Recline Leg Rest Elevation Rise		Bariatric Plus	
Basic Dimensions			
A Overall Height	46" / 1170mm		
B Back Height	29" / 730mm		
C Seat Height	20" / 510mm (4" Seat Cushion) 21" / 535mm (5" Seat Cushion)		
	18" / 450mm (4" Seat Cushion) 19" / 480mm (5" Seat Cushion)	Bariatric F	

Available Seat and Overall Widths

D Seat Width (SW)	E Overall Width	
26"/650mm	32"/820mm	
28"/700mm	34" / 870mm	
30"/750mm	36"/920mm	
34" / 850mm	45" / 1040mm	Bariatric Plus

Seat Depth

F 18 - 26" / 450mm - 650mm Adjustment achieved by removing fibre filling or zipping off lower waterfall cushion

Rise Angle

G 30° Forward Tilt & Rise

Seat Tilt Angle

Н	26° Backward Tilt	(650mm - 750mm SW on	ly)
Bariatric Plus does not have Tilt in Space Functionality		Bariatric	

Overall Length

I 48"/1200mm J 71"/1800mm

Patient Weight Limit

Ν	650mm - 750mm	40st / 254kg / 560lb	
Ν	850mm	70st / 450kg / 980lb	Ba

Leg Lift Capacity

	000	22007 110197 00010	Danamori
ı	850mm	22st / 140kg / 308lb	Bariatric Pl
	650mm - 750mm	20st / 12/kg / 280lb	

atric Plus

11

Chair Weight Unloaded

SW 650mm	12.8st / 81kg / 179lb	
SW 700mm	13st / 83kg / 183lb	
SW 750mm	13.4st / 85kg / 188lb	
SW 850mm	15st / 95kg / 210lb	Bariatric Plus

















Seating Matters are an **ISO 9001** Quality Management System approved company. © Seating Matters. Published March 2024



N. Ireland & ROI

131 Carnamuff Road Limavady Northern Ireland BT499JG T +44 (0) 28 777 666 24 contact@seatingmatters.com

Australia

9/256 New Line Road Dural NSW 2158 Australia T +611300001050 info@seatingmatters.com.au