Rotadorm Mobil

The transportable stand-up bed





- Mounted on the transport adapters, the bed is so small that it can be transported in almost any lift and is no longer held up by any door frame.
- **>** Four functions can be saved with the service handset.





The **Rotadorm Mobil** was developed for use in short-term care at home. The standup bed can be quickly and easily assembled and disassembled and stowed away to save space when it is not needed. Four individual settings can be saved using the service handset, which are a practical aid for patients and carers in everyday life.

Product details

- Mobile single bed with 90 cm width, adjustable to two lengths
- Under-bed clearance: 1 41 cm
- Four programmable positions: lying, turning, sitting and standing up
- 2-section special mattress with incontinence cover (optional)
- 5-section lying surface with wooden slats
- Patient lifting pole with grab handle (optional)
- Trendelenburg and anti-Trendelenburg positioning possible
- Additional stand-up function through programmable tilting of the thigh support
- Braking system on head and foot side
- Transport dimensions: $80 \times 96 \times 109 \text{ cm}$ (W x L x H)
- Transport adapter included
- Easy to transport and store thanks to small dimensions
- Removable hand rails

Aid numbers for German insurances: 19.40.01.8005 health ins. 50.45.01.4005 nursing care



Product images















The product images on this page show the standard version of the stand-up bed. The final product changes visually due to optional equipment ("Extra equipment/optional" on next page), such as an additional table. If you have any questions regarding an order, please do not hesitate to contact us \square .

pecial beds

Technical data



Mattress width: 90 cm

Mattress length: 200 / 210 cm

Bed height: 31 - 81 cm

Max. patient weight: 135 kg

Safe working load: 175 kg

Extra equipment / optional

Lying surface dimension: 90 x 200 cm (SI-013-1)

Lying surface dimension: 90 x 210 cm (SI-013-21)

Infrared control

Patient lifting pole with grab handle (NS-024-0)

provides secure and stable support to make it easier to sit up or change position in bed

Cover for handle receptacle (SI-011.80.725)

enables getting out of bed to the side

Table (SI-132-0)

serves as storage option during sitting posture, mounting on hand rail over wire frame, dimensions: 365 x 265 x 10 mm

Pillow clips (SI-038-0)

only needed for bed sheets in special colours – securing the pillow in position when the headrest is raised to aid comfort and convenience

Bed sheet - additional colours

see page 8 (pillow clips required)

Bedside cabinets (OM-011-B / OM-021-B)

bedside cabinets with table top, fixed height or height-adjustable


Alternating pressure system


SLK III digital - alternating pressure system


Mattress

Prime-Care (Standard)

Extra equipment / optional







Pillow clipsfor bed sheets in special colours


Table

Cover for handle receptacle

Mattress*

Available mattress






Description	Prime-Care	
Height	12 cm	
Pressure ulcer prophylaxis up to grade 2 (preventive)	√	
Material	SoliCore HRC	
Volume weight	45 kg/m³	
Visco-elastic	×	
Patient weight	up to 140 kg	

Alternating pressure system*

Alternating pressure system for decubitus prophylaxis and therapy

- Static and alternating pressure mode
- CPR function for resuscitation
- Large-cell mattress overlay system
- From low to high risk
- Optical low pressure alarm
- Patient weight: 40 195 kg
- Height: 21 cm



^{*} The size of the mattress and the alternating pressure system depends on the length of the respective bed. The right size is selected by the ISKOmed team.

Bedside cabinets*







Bedside cabinet with table top

(Dimensions: 34 x 57 x 88 cm) (Table top: 59 x 33 cm)







Bedside cabinet with table top (height-adjustable)

(Dimensions: approx. 34 x 68 x 88 cm) (Table top: 57 x 34 cm)

^{*}The colours of the bedside cabinets are perfectly matched to the Rotadorm Care and Mobil.

Bed sheets in special colours



(Pillow clips required)



Please note that due to technical limitations, the colours of the bed sheets cannot be displayed exactly in the swatches shown here. In addition, the use of different display technologies and your individual display settings can lead to distortions in the colour representation. The names of the colours serve an additional definition of the actual colour.



ISKOmed UK Limited · M +44 (0)7759 448 118 greg.whelan@iskomed.com · www.iskomed.com